




For additional help resources, scan this QR code to the left using your phone.

1. Open the Camera app from the Home screen or App menu.
2. Select the rear-facing camera and ensure that the code appears entirely in the view.
3. Tap the notification or pop-up to open the link associated with the QR code.
4. Alternatively, you can type myndimmersive.com/quickstart in any web browser to access resources.


STEP 1: Charge Headset

1. Plug charging cable into power adapter
2. Plug power adapter into wall outlet
3. Plug charging cable into headset; wait for green light (~2-2.5 hrs to full charge)




STEP 2: Power Devices

1. Unplug cable from headset
2. Hold power button for 3 sec
3. Remove battery activation tab from both controllers. Replace batteries (AA) as needed



STEP 3: Safety and Volume


1. Attach hand safety straps
2. Passthrough button
3. Volume button



STEP 4: Headset Fit

1. Put on headset
2. Adjust top strap
3. Use finger wheel to tighten headset

Note: Press and hold the infinity button (∞) on the controller for 2 sec to reorient the headset view.




STEP 5: Create Boundary

If prompted to create new boundary:


1. Ensure area is free of obstacles
2. Follow instructions in headset using trigger to make selections

Note: A boundary is a virtual play area that helps keep users safe within a clear space.



STEP 6: Finishing Session

1. Hold power button for 3 sec to power down headset. Light will turn off
2. Clean and sanitize
3. Plug charging cable into headset



Applications



Mynd Immersive 360 Video Player

Explore 200+ immersive videos designed to soothe, inspire, and engage.



Breathe XR

A guided breathing experience designed to support mindfulness and relaxation.



Progressive Muscle Relaxation (PMR)

A guided exercise that systematically contracts and relaxes muscle groups to help users recognize and reduce tension.



Soundscape Distraction

A customizable sound experience designed to help address discomfort and promote emotional balance.



Arrowhead - Medieval Archery VR

Step back in time to a medieval archery range, where your aim and precision are put to the ultimate test.



Colorway VR - Virtual Art Therapy

A recreational VR coloring experience that promotes creativity and focus.



Butterflies

Catch butterflies in a peaceful meadow using natural arm movements.



Puzzling Places

An immersive 3D puzzle experience that encourages focus, spatial reasoning, and problem-solving.



Power Pinball

A VR pinball experience that promotes visual tracking, timing, and responsive interaction in a dynamic format.



Bowling Blitz

A VR take on the beloved game of bowling that encourages movement, coordination, and focus through goal-driven play.